

PCR DIRECT® Guide to Product Selection

Arcis PCR DIRECT® Product	PCR DIRECT® Extracting Transport Medium	PCR DIRECT® Swab Extraction Reagent	PCR DIRECT® Saliva Extraction Reagent
For use at the point of collection*	✓		
Achieve greatest possible safety through earliest inactivation	✓		
For use in the Lab		✓	✓
Extract nucleic acids from dry samples. Avoid need for media		✓	✓
Simplify protocol: no incubation, no capture and concentrate, save time and money	✓	✓	✓
Reduce viscosity of samples, enable automation	✓	✓	✓
Achieve higher sensitivity compared to other direct methods	✓	✓	✓
Avoid guanidinium & hazardous waste	✓	✓	✓
Sequesters RNAses	✓	✓	✓
Lyophilisation compatible	✓	✓	✓
Sample Type	Swab/Saliva	Swab	Saliva

**Apply specified volume to sample with dropper or pre-aliquot into collection device.*